

抽吸导管的应用： **Why? When? How?**


颜红兵, MD, PhD, FACC, FESC

国家心血管病中心 中国医学科学院阜外医院

抽吸导管的应用

■ Why?


Proximal



Guiding Catheter:

- inexpensive
- less selective
- systemic effects

Intralesional



ClearWay Catheter:

- additional costs
- controlled and targeted drug delivery

Distal



Microcatheters:

- thrombectomy devices
- high availability
- Infusion in microcirculatory bed

抽吸导管：指南建议

Title	Citation		Class	LOE
2012 ESC Guidelines ST-segment elevation myocardial infarction . 	European Heart Journal 2012 Oct;33(20):2569-619	Routine aspiration should be considered	IIa	B
2014 ESC/EACTS guidelines on myocardial revascularization 	Eur Heart J. 2014 Oct 1;35(37):2541-619	May be considered in selected patients	IIb	A
2015 ACC/AHA focused update PPCI 	JACC	Routine thrombectomy not useful	III	A
2015 ACC/AHA focused update PPCI 	JACC	Selective and bailout Thrombectomy not well established	IIb	C
2017 ESC Guidelines ST-segment elevation myocardial infarction 	European Heart Journal 2017	Routine use of thrombus aspiration is not recommended.	III	A

单纯抽吸就能降低死亡率吗？



What is new in 2017 Guidelines on AMI-STEMI




2012	CHANGE IN RECOMMENDATIONS	2017
	Radial access	MATRIX
	DES over BMS	EXAMINATION, COMFORTABLE-AMI, NORSTENT
	Complete Revascularisation	PRAMI, DANAMI-3-PRIMULTI, CVLPRIT, Compare-Acute
	Thrombus Aspiration	TOTAL, TASTE
	Bivalirudin	MATRIX, HEAT-PPCI
	Enoxaparin	ATOLL, Meta-analysis
	Early Hospital Discharge	Small trials & observational data
Oxygen when SaO ₂ <95%	OXYGEN	Oxygen when SaO ₂ <90% AVOID, DETO2X
Same dose i.v. in all patients	TNK-tPA	Half dose i.v. in Pts ≥75 years STREAM

“3T”研究说明了什么？

- TAPAS (2008)
- TASTE (2013)
- TOTAL (2015)



指南与实践




面临的主要挑战

- Shorten time to reperfusion
- Reduce access site related complications
- Choosing the correct stent
- High thrombotic burden
- No reflow phenomenon
- Reperfusion injury
- STEMI with multivessel disease
- STEMI with cardiogenic shock
- Long-term DAPT management




无复流发生机制

1. Distal embolization



2. Ischemia-related injury 3. Reperfusion-related injury



4. Individual susceptibility

- Genetic variability
- Diabetes
- Acute Hyperglycemia
- Dislipidemia
- Lack of preconditioning


adpted from Niccoli et al. J Am Coll Cardiol 2009;54:281-92

ACS患者的特征

- **高炎症状态**
- **高致血栓状态**和血栓成分的复杂性
- 合并多支病变率高（60.9%）
- 合并性疾病率高（糖尿病、肾功能不全和房颤等）

Lessons From Sudden Coronary Death A Comprehensive Morphological Classification Scheme for Atherosclerotic Lesions


Renu Virmani, Frank D. Kolodgie, Allen P. Burke, Andrew Farb, Stephen M. Schwartz



这3种类型对临床影响一样吗？

■ 要点2:

同样的斑块类型风险不一样，还取决于斑块下的化学结构！




■ 要点1:

斑块病理类型不同对临床的影响可能不同，似乎斑块破裂的风险更高





■ 要点3:

同样斑块下结构风险不一样，还取决于炎症水平



抽吸物的组成

- 血小板、新鲜血栓、机化的组织斑块
- 动脉粥样硬化斑块组成
- 胆固醇结晶、泡沫细胞、坏死的核心
- 血管活性物质
- 5羟色胺、ADP、血栓素A₂
- 炎性介质




“血栓”的动态性质

Culprit Vessel

Fresh, Friable Platelet
rich Thrombus

Dense Organized
Thrombus with RBC,
Fibrin





- Aspirin
- Ticagrelor
- Clopidogrel
- Bivalirudin
- Abciximab

- Thrombo – Aspiration
- Angiojet (Large Vessel)
- Stenting

- Angiojet
- Clearway
- Laser
- Staged procedure with prolonged anti-coagulation


AMI
presentation
in developing
countries

STEMI患者血栓分层处理策略

分级	血栓定义	血管造影举例	Mehta分类	技术要点	
0	No cine angiographic characteristics of thrombus present			Aspiration Catheter	Angio Jet
1	Possible thrombus present. Angiography demonstrates reduced contrast density, haziness, irregular lesion contour or a smooth convex "meniscus" at the site of total occlusion suggestive but not diagnostic of thrombus		Direct Stent +/- Pre-dilatation	•Most effective with fresh clot; organized thrombus is more resistant to debulking. •Have different profiles, different push-ability, tractability and aspiration rates. •All are 6F-compatible. It is useful to stock and be familiar with the use of at least one. •Flush catheter lumen well before use as it facilitates better tracking over the wire.	•Can be used from the radial route. Although LAD and some LCX may not need a TPM, I place TPM's in all Angiojet procedures. •Often, multiple passes will be required. Try to pause after every 2-3 passes to enable hemodynamics to be restored, to optimize guide wire and guiding catheter support and to evaluate the results.
2	Thrombus present-small size: Definite thrombus with greatest dimensions less than or equal to 1/2 vessel diameter			Aspiration thrombectomy	•Often, just the first passage will restore adequate flow •Resistant and stubborn thrombus will require more distal advancement that must be done more carefully.
3	Thrombus present- moderate size: Definite thrombus but with greatest linear dimension greater than 1/2 but less than 2 vessel diameters			•Avoid kinking the catheter – advance slowly over the initial, softer portion of the catheter.	•Avoid advancing in severe tortuosity and in vessels < 2mm
4	Thrombus present- large size: As in Grade 3 but with the largest dimension greater than or equal to 2 vessel diameters			Angio Jet	•Monitor distal tip of the guide wire as the aspiration catheter is advanced – it is not uncommon for the guide wire to advance during this maneuver
5	Total occlusion			•Advance the aspiration catheter through the entire length of occlusive disease.	•Since the Angiojet is used for large thrombus burden and high thrombus grade, consider Abciximab as adjunctive therapy

为何使用抽吸导管？

- Restore flow in the culprit artery.
- Optimize myocardial perfusion.




抽吸导管只是一个平台

- 可以抽吸
- 可以超选择给药
- 可以超选择造影
- 各种性质的血栓
- 斑块物质
- 降低局部炎性介质

抽吸导管的应用

■ When?


Proximal



Guiding Catheter:

- inexpensive
- less selective
- systemic effects

Intralesional



ClearWay Catheter:

- additional costs
- controlled and targeted drug delivery

Distal




Microcatheters:

- thrombectomy devices
- high availability
- Infusion in microcirculatory bed

你的目的？


- 抽吸？
- 超选择给药？
- 超选择造影？



抽吸导管的应用

■ How?


Proximal



Guiding Catheter:

- inexpensive
- less selective
- systemic effects

Intralesional



ClearWay Catheter:

- additional costs
- controlled and targeted drug delivery

Distal




Microcatheters:

- thrombectomy devices
- high availability
- Infusion in microcirculatory bed

抽吸导管的应用

✓ 抽吸导管具有3大功能：

- ✓ 抽吸
- ✓ 输送
- ✓ 支撑




✓ 抽吸导管只是提供了一个平台

✓ 终点是要改善血流与灌注

✓ 应当强调综合处理

抽吸导管的应用

✓ 请提问题



抽吸导管的应用

- 不忘初心

改善血流与灌注

- 牢记使命

各种性质的血栓

斑块物质

降低局部炎性介质